



**TITAN**  
REGULAR TEMPER  
MOLAR BANDS

**TITAN**  
FIRM TEMPER  
MOLAR BANDS

**MAXUM**  
ANATOMICAL  
MOLAR BANDS

**SIERRA**  
GENERAL PURPOSE  
MOLAR BANDS

**AERO**  
*Precision Fit*  
*Bicuspid Bands*

**FAIRFIELD**  
Prewelded Bands

## A note on our Bands...

Fairfield Orthodontics offers a wide range of first and second molar bands, bicuspid bands, and cuspid and anterior bands. We have regular temper and firm temper, anatomically fitted, as well as wide and narrow molar bands to suit the doctor's every need.

Our popular Titan Regular Temper molar bands are available in half sizes 29.5 through 44. Sizes are permanently marked for easy identification. They are anatomically contoured for an accurate fit and like all of our bands are available etched for better retention.

Fairfield Orthodontics is unique in that we inventory prewelded molar bands. If an orthodontic practice places their trust in us as their supplier for molar or bicuspid bands, we are committed to provide them with the fastest service possible. We inventory every doctor's prescription, so it is always in stock and ready to ship.

Fairfield Orthodontics wants to make the ordering process as simple as possible. Therefore, we developed the ez fax ordering system. On the back of this page is a fax order form which can be copied and used for placing band orders. You will find a similar order form at the beginning of each section throughout the catalog- each specifically designed for ordering the products in that section.

**Bill To:**

NAME: \_\_\_\_\_

ADDRESS 1: \_\_\_\_\_

ADDRESS 2: \_\_\_\_\_

CITY: \_\_\_\_\_ STATE: \_\_\_\_\_ ZIP: \_\_\_\_\_

PHONE: ( ) \_\_\_\_\_

FAX: ( ) \_\_\_\_\_

**Ship To:**  Same as Bill To

NAME: \_\_\_\_\_

ADDRESS 1: \_\_\_\_\_

ADDRESS 2: \_\_\_\_\_

CITY: \_\_\_\_\_ STATE: \_\_\_\_\_ ZIP: \_\_\_\_\_

PHONE: ( ) \_\_\_\_\_

FAX: ( ) \_\_\_\_\_

**Band Form 1:**

Plain Bands     1st Molar     Etched

Prewelded Bands     2nd Molar

Upper Buccal Tube: \_\_\_\_\_

Upper Lingual Attach: \_\_\_\_\_

Lower Buccal Tube: \_\_\_\_\_

Lower Lingual Attach: \_\_\_\_\_

Special Instructions: \_\_\_\_\_

**Other Items:**



Qty	Item #	Description

**Shipping Options:**

 DHL Ground Service     Other \_\_\_\_\_

**Payment Options:**

Bill Me Net 30 Days (Overdue Fees May Apply)

Credit Card Payment     

\_\_\_\_\_|\_\_\_\_\_|\_\_\_\_\_|\_\_\_\_\_|\_\_\_\_\_|\_\_\_\_\_|\_\_\_\_\_|\_\_\_\_\_|\_\_\_\_\_|\_\_\_\_\_|

Sign. \_\_\_\_\_ Exp. \_\_\_\_/\_\_\_\_

**Band Form 2:**

Maxum	TITAN	Upper		Lower	
		Right	Left	Right	Left
1	29.5				
2	30				
3	30.5				
4	31				
5	31.5				
6	32				
7	32.5				
8	33				
9	33.5				
10	34				
11	34.5				
12	35				
13	35.5				
14	36				
15	36.5				
16	37				
17	37.5				
18	38				
19	38.5				
20	39				
21	39.5				
22	40				
23	40.5				
24	41				
25	41.5				
26	42				
27	42				
28	42				
29	42.5				
30	42.5				
31	43				
32	43				
33	43.5				
34	43.5				
35	44				
36	44				



# Bands

## TITAN REGULAR TEMPER MOLAR BANDS

Titan molar bands are designed to seat out of occlusion and are fully proportioned in all dimensions - diameter, height, and thickness - to provide maximum strength and fit with minimum trimming.

These regular temper bands are made more maleable by a heat treat process, allowing close anatomical adaptation.

Broad buccal and lingual surfaces are provided for superior retention and allow flexibility in positioning attachments more towards the occlusal or gingival, as desired.



**Micro-Etching**  
available on all bands

**Anatomical Design** for fast placement and accurate fit

**Straight Interproximal Wall**  
makes band seating easy



**Permanent Marking** for easy identification.

**Smooth Edges** for optimal patient comfort

**Ideally Tempered** for easy contouring

## TITAN REGULAR TEMPER MOLAR BANDS Maxillary 1<sup>st</sup> & 2<sup>nd</sup> Molar Band Kits

	29.5	30	30.5	31	31.5	32	32.5	33	33.5	34	34.5	35	35.5	36	36.5	37	37.5	38	38.5	39	39.5	40	40.5	41	41.5	42	42.5	43	43.5	44
<b>Trial Kit</b> (80 molar bands)	-	-	-	-	-	2	2	4	4	4	6	6	10	10	10	6	6	4	4	2	-	-	-	-	-	-	-	-	-	-
<b>Intro Kit</b> (150 molar bands)	2	2	2	2	2	4	4	6	6	8	10	10	10	12	12	10	8	8	6	4	4	2	2	2	2	2	2	2	2	2
<b>Professional Kit</b> (300 molar bands)	4	4	4	4	6	8	8	10	12	16	16	20	24	24	24	20	16	16	12	10	8	6	6	6	4	4	2	2	2	2

See our  
prewelded  
bands on  
page 45

**Trial Kit** (80 molar bands)

#144-500 Maxillary First  
#144-520 Maxillary Second

**Intro Kit** (150 molar bands)

#144-600 Maxillary First  
#144-620 Maxillary Second

**Pro Kit** (300 molar bands)

#144-700 Maxillary First  
#144-720 Maxillary Second

All kits contain half right and half left. Distribution may be altered as long as total number of bands remains the same.

## TITAN Mandibular 1<sup>st</sup> Molar Band Kits

REGULAR TEMPER  
MOLAR BANDS

	29.5	30	30.5	31	31.5	32	32.5	33	33.5	34	34.5	35	35.5	36	36.5	37	37.5	38	38.5	39	39.5	40	40.5	41	41.5	42	42.5	43	43.5	44
<b>Trial Kit</b> (80 molar bands)	-	-	-	-	-	2	2	4	4	4	6	6	10	10	10	6	6	4	4	2	-	-	-	-	-	-	-	-	-	-
<b>Intro Kit</b> (150 molar bands)	2	2	2	2	2	4	4	6	6	8	10	10	10	12	12	10	8	8	6	4	4	2	2	2	2	2	2	2	2	2
<b>Professional Kit</b> (300 molar bands)	4	4	4	4	6	8	8	10	12	16	16	20	24	24	24	20	16	16	12	10	8	6	6	6	4	4	2	2	2	2

See our  
prewelded  
bands on  
page 45

**Trial Kit** (80 molar bands)  
#144-550 Mandibular First

**Intro Kit** (150 molar bands)  
#144-650 Mandibular First

**Pro Kit** (300 molar bands)  
#144-750 Mandibular First

All kits contain half right and half left. Distribution may be altered as long as total number of bands remains the same.

## TITAN Mandibular 2<sup>nd</sup> Molar Band Kits

REGULAR TEMPER  
MOLAR BANDS

	29.5	30	30.5	31	31.5	32	32.5	33	33.5	34	34.5	35	35.5	36	36.5	37	37.5	38	38.5	39	39.5	40	40.5	41	41.5	42	42.5	43	43.5	44
<b>Trial Kit</b> (80 molar bands)	-	-	-	-	-	2	2	4	4	4	6	6	10	10	10	6	6	4	4	2	-	-	-	-	-	-	-	-	-	-
<b>Intro Kit</b> (150 molar bands)	2	2	2	2	2	4	4	6	6	8	10	10	10	12	12	10	8	8	6	4	4	2	2	2	2	2	2	2	2	2
<b>Professional Kit</b> (300 molar bands)	4	4	4	4	6	8	8	10	12	16	16	20	24	24	24	20	16	16	12	10	8	6	6	6	4	4	2	2	2	2

See our  
prewelded  
bands on  
page 45

**Trial Kit** (80 molar bands)  
#144-570 Mandibular Second

**Intro Kit** (150 molar bands)  
#144-670 Mandibular Second

**Pro Kit** (300 molar bands)  
#144-770 Mandibular Second

All kits contain half right and half left. Distribution may be altered as long as total number of bands remains the same.

TITAN  
REGULAR TEMPER  
MOLAR BANDS



We stock all our prewelded bands so they are always ready to ship the same day!

Details on page 57

# Bands

## TITAN FIRM TEMPER MOLAR BANDS

Titan molar bands are designed to seat out of occlusion and are fully proportioned in all dimensions - diameter, height, and thickness - to provide maximum strength and fit with minimum trimming.

These firm temper bands provide a "spring fit" for excellent mechanical retention.

Broad buccal and lingual surfaces are provided for superior retention and allow flexibility in positioning attachments more towards the occlusal or gingival, as desired.



**Micro-Etching**  
available on all bands

**Anatomical Design** for fast placement and accurate fit

**Straight Interproximal Wall**  
makes band seating easy



**Permanent Marking** for easy identification.

**Smooth Edges** for optimal patient comfort

**Firm Temper** for spring fit and excellent retention

## TITAN FIRM TEMPER MOLAR BANDS Maxillary 1<sup>st</sup> & 2<sup>nd</sup> Molar Band Kits

	129.5	130	130.5	131	131.5	132	132.5	133	133.5	134	134.5	135	135.5	136	136.5	137	137.5	138	138.5	139	139.5	140	140.5	141	141.5	142	142.5	143	143.5	144
<b>Trial Kit</b> (80 molar bands)	-	-	-	-	-	2	2	4	4	4	6	6	10	10	10	6	6	4	4	2	-	-	-	-	-	-	-	-	-	-
<b>Intro Kit</b> (150 molar bands)	2	2	2	2	2	4	4	6	6	8	10	10	10	12	12	10	8	8	6	4	4	2	2	2	2	2	2	2	2	2
<b>Professional Kit</b> (300 molar bands)	4	4	4	4	6	8	8	10	12	16	16	20	24	24	24	20	16	16	12	10	8	6	6	6	4	4	2	2	2	2

**FIRM TEMPER**

**Trial Kit** (80 molar bands)

#140-500 Maxillary First  
#140-520 Maxillary Second

**FIRM TEMPER**

**Intro Kit** (150 molar bands)

#140-600 Maxillary First  
#140-620 Maxillary Second

**FIRM TEMPER**

**Pro Kit** (300 molar bands)

#140-700 Maxillary First  
#140-720 Maxillary Second

See our  
prewelded  
bands on  
page 45

All kits contain half right and half left. Distribution may be altered as long as total number of bands remains the same.

## TITAN Mandibular 1<sup>st</sup> Molar Band Kits

	129.5	130	130.5	131	131.5	132	132.5	133	133.5	134	134.5	135	135.5	136	136.5	137	137.5	138	138.5	139	139.5	140	140.5	141	141.5	142	142.5	143	143.5	144
<b>Trial Kit</b> (80 molar bands)	-	-	-	-	-	2	2	4	4	4	6	6	10	10	10	6	6	4	4	2	-	-	-	-	-	-	-	-	-	-
<b>Intro Kit</b> (150 molar bands)	2	2	2	2	2	4	4	6	6	8	10	10	10	12	12	10	8	8	6	4	4	2	2	2	2	2	2	2	2	2
<b>Professional Kit</b> (300 molar bands)	4	4	4	4	6	8	8	10	12	16	16	20	24	24	24	20	16	16	12	10	8	6	6	6	4	4	2	2	2	2

See our  
prewelded  
bands on  
page 45

### FIRM TEMPER

**Trial Kit** (80 molar bands)  
#140-550 Mandibular First

### FIRM TEMPER

**Intro Kit** (150 molar bands)  
#140-650 Mandibular First

### FIRM TEMPER

**Pro Kit** (300 molar bands)  
#140-750 Mandibular First

All kits contain half right and half left. Distribution may be altered as long as total number of bands remains the same.

## TITAN Mandibular 2<sup>nd</sup> Molar Band Kits

Note: Titan Lower Second Molar Bands are only available Regular Temper.

	29.5	30	30.5	31	31.5	32	32.5	33	33.5	34	34.5	35	35.5	36	36.5	37	37.5	38	38.5	39	39.5	40	40.5	41	41.5	42	42.5	43	43.5	44
<b>Trial Kit</b> (80 molar bands)	-	-	-	-	-	2	2	4	4	4	6	6	10	10	10	6	6	4	4	2	-	-	-	-	-	-	-	-	-	-
<b>Intro Kit</b> (150 molar bands)	2	2	2	2	2	4	4	6	6	8	10	10	10	12	12	10	8	8	6	4	4	2	2	2	2	2	2	2	2	2
<b>Professional Kit</b> (300 molar bands)	4	4	4	4	6	8	8	10	12	16	16	20	24	24	24	20	16	16	12	10	8	6	6	6	4	4	2	2	2	2

See our  
prewelded  
bands on  
page 45

**Trial Kit** (80 molar bands)  
#144-570 Mandibular Second

**Intro Kit** (150 molar bands)  
#144-670 Mandibular Second

**Pro Kit** (300 molar bands)  
#144-770 Mandibular Second

All kits contain half right and half left. Distribution may be altered as long as total number of bands remains the same.

**TITAN**  
FIRM TEMPER  
MOLAR BANDS



We stock all our prewelded bands so they are always ready to ship the same day!

Details on page 57

# Bands

## MAXUM ANATOMICAL MOLAR BANDS

Maxum Anatomical molar bands are designed to provide a "snap on" fit that provides excellent mechanical retention and requires minimum custom contouring.

Consistent graduation of sizes and specific right and left upper and lower anatomy eliminates the need for cutting and lapping.

The contoured design eliminates gingival overhang. Maxum also features a hardened occlusal perimeter that resists deformation during band seating.

Available in Maxillary First, Mandibular First, and Mandibular Second Anatomies.



**Micro-Etching**  
available on all bands

**Anatomical Design** for fast placement and accurate fit

**Hardened occlusal perimeter** resists deformation during seating



**Permanent Marking** for easy identification.

**Smooth Edges** for optimal patient comfort

**Ideally Tempered** for easy contouring

### MAXUM ANATOMICAL MOLAR BANDS

## Maxillary 1<sup>st</sup> & 2<sup>nd</sup> Molar Band Kits

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36
<b>Trial Kit</b> (80 molar bands)	-	-	-	-	-	-	-	-	-	-	-	-	2	2	2	2	2	4	6	6	6	6	6	6	6	6	6	6	4	2	-	-	-	-	-	-
<b>Intro Kit</b> (150 molar bands)	2	2	2	2	2	2	2	2	2	2	2	2	4	4	4	4	4	4	8	8	8	8	8	8	8	8	8	8	4	4	4	2	2	2	2	2
<b>Professional Kit</b> (300 molar bands)	2	2	2	4	4	4	4	4	4	6	6	6	8	8	8	8	10	12	14	16	16	16	16	16	16	16	14	12	10	8	6	6	4	4	4	4

See our  
prewelded  
bands on  
page 45

**Trial Kit** (80 molar bands)

#150-500 Maxillary First  
#150-520 Maxillary Second

**Intro Kit** (150 molar bands)

#150-600 Maxillary First  
#150-620 Maxillary Second

**Pro Kit** (300 molar bands)

#150-700 Maxillary First  
#150-720 Maxillary Second

All kits contain half right and half left. Distribution may be altered as long as total number of bands remains the same.





## Mandibular 1<sup>st</sup> Molar Band Kits

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36
<b>Trial Kit</b> (80 molar bands)	-	-	-	-	-	-	-	-	-	-	-	-	2	2	2	2	2	4	6	6	6	6	6	6	6	6	6	6	4	2	-	-	-	-	-	-
<b>Intro Kit</b> (150 molar bands)	2	2	2	2	2	2	2	2	2	2	2	2	4	4	4	4	4	4	8	8	8	8	8	8	8	8	8	8	4	4	4	2	2	2	2	2
<b>Professional Kit</b> (300 molar bands)	2	2	2	4	4	4	4	4	4	6	6	6	8	8	8	8	10	12	14	16	16	16	16	16	16	16	14	12	10	8	6	6	4	4	4	4

See our  
prewelded  
bands on  
page 45

**Trial Kit** (80 molar bands)  
#150-550 Mandibular First

**Intro Kit** (150 molar bands)  
#150-650 Mandibular First

**Pro Kit** (300 molar bands)  
#150-750 Mandibular First

All kits contain half right and half left. Distribution may be altered as long as total number of bands remains the same.



## Mandibular 2<sup>nd</sup> Molar Band Kits

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32		
<b>Trial Kit</b> (80 molar bands)	-	-	-	-	-	-	-	-	2	2	2	4	4	4	6	6	6	6	6	6	6	6	6	6	2	-	-	-	-	-	-	-	-	
<b>Intro Kit</b> (150 molar bands)	2	2	2	2	2	2	2	2	4	4	4	4	4	4	4	10	10	10	10	10	10	10	10	10	2	2	2	2	2	2	2	2	2	2
<b>Professional Kit</b> (300 molar bands)	4	4	4	4	4	4	4	4	8	8	8	8	8	8	8	20	20	20	20	20	20	20	20	20	4	4	4	4	4	4	4	4	4	4


See our  
prewelded  
bands on  
page 45

**Trial Kit** (80 molar bands)  
#150-570 Mandibular Second

**Intro Kit** (150 molar bands)  
#150-670 Mandibular Second

**Pro Kit** (300 molar bands)  
#150-770 Mandiublar Second

All kits contain half right and half left. Distribution may be altered as long as total number of bands remains the same.



We stock all our prewelded bands so they are always ready to ship the same day!

Details on page 57

# Bands

## SIERRA GENERAL PURPOSE MOLAR BANDS

Sierra molar bands offer an exceptionally good fit on upper or lower, first or second molars. Sizes are proportionally graduated, with bold permanent marking for easy identification.

Sierra molar bands allow for reduced inventory due to composite right/left anatomy yet offer many features. Sierra bands are optimal temper which adapt easily to a variety of tooth shapes, feature a buccal groove for seating reference and are available in WIDE (.25") or NARROW (.150").



Available Maxillary Wide or Narrow and Mandibular Wide or Narrow. (Narrow pictured above)

**Micro-Etching**  
available on all bands

**Permanent Marking** for  
easy identification.

**Smooth Edges** for  
optimal patient comfort



**SIERRA**

**Composite left/right anatomy**  
reduces inventory

**Buccal Groove** provides  
seating reference

**Ideally Tempered** for  
easy contouring

## SIERRA Wide Maxillary 1<sup>st</sup> & 2<sup>nd</sup> Molar Band Kits

GENERAL PURPOSE  
MOLAR BANDS

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36
<b>Trial Kit</b> (80 molar bands)	-	-	-	-	-	-	-	-	-	-	-	-	2	2	2	2	2	4	6	6	6	6	6	6	6	6	6	6	4	2	-	-	-	-	-	-
<b>Intro Kit</b> (150 molar bands)	2	2	2	2	2	2	2	2	2	2	2	2	4	4	4	4	4	4	8	8	8	8	8	8	8	8	8	8	4	4	4	2	2	2	2	2
<b>Professional Kit</b> (300 molar bands)	2	2	2	4	4	4	4	4	4	6	6	6	8	8	8	8	10	12	14	16	16	16	16	16	16	16	16	14	12	10	8	6	6	4	4	4

See our  
prewelded  
bands on  
page 45

**Trial Kit** (80 molar bands)  
#160-500 Maxillary Wide

**Intro Kit** (150 molar bands)  
#160-600 Maxillary Wide

**Pro Kit** (300 molar bands)  
#160-700 Maxillary Wide

All kits contain half right and half left. Distribution may be altered as long as total number of bands remains the same.

## SIERRA *Wide* Mandibular 1<sup>st</sup> & 2<sup>nd</sup> Molar Band Kits

GENERAL PURPOSE  
MOLAR BANDS

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36
<b>Trial Kit</b> (80 molar bands)	-	-	-	-	-	-	-	-	-	-	-	-	2	2	2	2	2	4	6	6	6	6	6	6	6	6	6	6	4	2	-	-	-	-	-	-
<b>Intro Kit</b> (150 molar bands)	2	2	2	2	2	2	2	2	2	2	2	2	4	4	4	4	4	4	8	8	8	8	8	8	8	8	8	8	4	4	4	2	2	2	2	2
<b>Professional Kit</b> (300 molar bands)	2	2	2	4	4	4	4	4	4	4	6	6	6	8	8	8	8	10	12	14	16	16	16	16	16	16	16	14	12	10	8	6	6	4	4	4

See our  
prewelded  
bands on  
page 45

**Trial Kit** (80 molar bands)  
#160-550 Mandibular Wide

**Intro Kit** (150 molar bands)  
#160-650 Mandibular Wide

**Pro Kit** (300 molar bands)  
#160-750 Mandibular Wide

## SIERRA *Narrow* Maxillary 1<sup>st</sup> & 2<sup>nd</sup> Molar Band Kits

GENERAL PURPOSE  
MOLAR BANDS

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36
<b>Trial Kit</b> (80 molar bands)	-	-	-	-	-	-	-	-	-	-	-	-	2	2	2	2	2	4	6	6	6	6	6	6	6	6	6	6	4	2	-	-	-	-	-	-
<b>Intro Kit</b> (150 molar bands)	2	2	2	2	2	2	2	2	2	2	2	2	4	4	4	4	4	4	8	8	8	8	8	8	8	8	8	8	4	4	4	2	2	2	2	2
<b>Professional Kit</b> (300 molar bands)	2	2	2	4	4	4	4	4	4	4	6	6	6	8	8	8	8	10	12	14	16	16	16	16	16	16	16	14	12	10	8	6	6	4	4	4

See our  
prewelded  
bands on  
page 45

**Trial Kit** (80 molar bands)  
#163-500 Maxillary Narrow

**Intro Kit** (150 molar bands)  
#163-600 Maxillary Narrow

**Pro Kit** (300 molar bands)  
#163-700 Maxillary Narrow

## SIERRA *Narrow* Mandibular 1<sup>st</sup> & 2<sup>nd</sup> Molar Band Kits

GENERAL PURPOSE  
MOLAR BANDS

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36
<b>Trial Kit</b> (80 molar bands)	-	-	-	-	-	-	-	-	-	-	-	-	2	2	2	2	2	4	6	6	6	6	6	6	6	6	6	6	4	2	-	-	-	-	-	-
<b>Intro Kit</b> (150 molar bands)	2	2	2	2	2	2	2	2	2	2	2	2	4	4	4	4	4	4	8	8	8	8	8	8	8	8	8	8	4	4	4	2	2	2	2	2
<b>Professional Kit</b> (300 molar bands)	2	2	2	4	4	4	4	4	4	4	6	6	6	8	8	8	8	10	12	14	16	16	16	16	16	16	16	14	12	10	8	6	6	4	4	4

See our  
prewelded  
bands on  
page 45

**Trial Kit** (80 molar bands)  
#163-550 Mandibular Narrow

**Intro Kit** (150 molar bands)  
#163-650 Mandibular Narrow

**Pro Kit** (300 molar bands)  
#163-750 Mandibular Narrow

# Bands

## **AERO** Precision Fit Bicuspid Bands

AERO precision fit bicuspid bands are known for their exceptional fit and retention. Three distinctive band shapes, and ideal contour assure proper fit and prevent the band from seating below the height of the contour.

Available in a wide selection of proportionally graduated sizes in specific anatomy for upper bicuspid (1st and 2nd), lower 1st bicuspid, and lower 2nd bicuspid.



**Lower 1st Bicuspid**

**Upper 1st & 2nd Bicuspid**

**Lower 2nd Bicuspid**



\*Many prewelded bicuspid band options available.  
Please call for details.

## **AERO** Precision Fit Bicuspid Bands

### Maxillary 1<sup>st</sup> & 2<sup>nd</sup> Bicuspid

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26
<b>Trial Kit</b> (80 bicuspid bands)	-	-	-	-	-	-	4	4	8	12	12	12	12	8	4	4	-	-	-	-	-	-	-	-	-	-
<b>Intro Kit</b> (200 bicuspid bands)	2	2	2	2	10	10	10	10	16	16	16	16	16	16	10	10	10	10	6	4	4	2	2	2	-	-
<b>Professional Kit</b> (350 bicuspid bands)	4	4	4	4	16	16	20	20	30	30	30	30	30	20	18	18	18	14	4	4	4	4	4	4	-	-

**Trial Kit** (80 bands)

#134-500 Maxillary Bicuspid

**Intro Kit** (200 bands)

#134-600 Maxillary Bicuspid

**Pro Kit** (350 bands)

#134-700 Maxillary Bicuspid

All kits contain half right and half left. Distribution may be altered as long as total number of bands remains the same.

# Bands

**AERO**  
Precision Fit  
Bicuspid Bands

## Mandibular 1<sup>st</sup> Bicuspid

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26
<b>Trial Kit</b> (40 bicuspid bands)	-	-	-	-	-	-	2	2	4	6	6	6	6	4	2	2	-	-	-	-	-	-	-	-	-	-
<b>Intro Kit</b> (78 bicuspid bands)	2	2	2	2	2	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4	2	2	2	2	-	-
<b>Professional Kit</b> (130 bicuspid bands)	4	4	4	4	6	6	6	6	6	6	6	6	6	6	6	6	6	6	6	6	6	4	4	4	-	-

**Trial Kit** (40 bands)  
#134-550 Mandibular First

**Intro Kit** (78 bands)  
#134-650 Mandibular First

**Pro Kit** (130 bands)  
#134-750 Mandibular First

All kits contain half right and half left. Distribution may be altered as long as total number of bands remains the same.

**AERO**  
Precision Fit  
Bicuspid Bands

## Mandibular 2<sup>nd</sup> Bicuspid

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26
<b>Trial Kit</b> (60 bicuspid bands)	-	-	-	-	-	-	-	-	2	4	4	10	10	10	10	4	4	2	-	-	-	-	-	-	-	-
<b>Intro Kit</b> (122 bicuspid bands)	-	-	2	2	2	2	4	6	6	6	8	8	10	10	8	8	6	6	6	6	4	4	2	2	2	2
<b>Professional Kit</b> (220 bicuspid bands)	-	-	4	4	4	6	10	10	10	10	14	14	18	18	14	14	10	10	10	10	6	6	6	4	4	4

**Trial Kit** (60 bands)  
#134-500 Mandibular Second

**Intro Kit** (122 bands)  
#134-600 Mandibular Second

**Pro Kit** (220 bands)  
#134-700 Mandibular Second

All kits contain half right and half left. Distribution may be altered as long as total number of bands remains the same.

**Prewelded Bicuspid  
Bands Available**



We have a large selection of bicuspid bands with attachments prewelded IN STOCK and READY SHIP. Please call for details.

**Also Available  
Anterior & Cuspid Bands**



